

β -Tubulin Mouse mAb Catalog NO.: EM1029 For research use only.

Overview

Product name β -Tubulin Mouse Monoclonal antibody

Source Mouse

Applications WB IHC IF

Species reactivity Human Rat Mouse Monkey Dog Chicken Hamster RabbitSheep Insect

Yeast

Recommended dilutions WesternBlot:1/5000

Immunohistochemistry:1/200 Immunofluorescence:1/100-200

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species Human

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band 55kDa

GenelD (Human) 203068

Human Swiss-Prot No. P07437

Cellular localization CytoplasmCytoskeletonMicrotubule


Alternative Names N/A

Background Microtubules are constituent parts of the mitotic apparatus cilia flagella and


elements of the cytoskeleton. They consist principally of 2 soluble proteins alpha- and beta-tubulin each of about 55000 Da. Antibodies against beta Tubulin are useful as loading controls for Western Blotting. However it should be noted that levels of $\,\beta$ -Tubulin may not be stable in certain cells.

https://www.elkbiotech.com T:86-27-59760950 E: elkbio@elkbiotech.com


For example expression of β -Tubulin in adipose tissue is very low and therefore β -Tubulin should not be used as loading control for these tissues.




Western blot analysis of A549 (1) Rat brain (2) Mouse brain (3) Chicken lung (4) and Rabbit testis (5)Sheep muscle (6) with β -tubulin mouse mAb(5G3) diluted at:5000.



IHC Staining of Human colon tissue with β-tubulin mouse mAb(5G3) diluted at:200.



IF analysis of Hela with EM1029(Left) and DAPI (Right) diluted at:100.



Western blot analysis of extracts from HeLa (Lane)MCF-7 (Lane 2) MG63 (Lane 3) Min6 (Lane 4)Mouse Kindey (Lane 5) Rat Kindey (Lane 6)Rat Brain (Lane 7) using β -tubulin diluted at 1:10000.