

NIF3L rabbit pAb

Cat No.:ES14492

For research use only

Overview

Product Name	NIF3L rabbit pAb
Host species	Rabbit
-	
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human NIF3L AA
	range: 85-135
Specificity	This antibody detects endogenous levels of NIF3L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	NIF3L
Gene Name	NIF3L1 ALS2CR1 MDS015 My018
Cellular localization	Cytoplasm . Nucleus . Interaction with COPS2 may
	regulate localization to the nucleus
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	60491
Human Swiss-Prot Number	Q9GZT8
Alternative Names	
Background	




ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C