

NIBAN rabbit pAb

Cat No.:ES14493

For research use only

Overview

Product Name	NIBAN rabbit pAb	
Host species	Rabbit	
Applications	WB	
Species Cross-Reactivity	Human; Mouse;Rat	
Recommended dilutions	WB 1: 500-2000	
Immunogen	Synthesized peptide derived from human NIBAN AA range: 334-384	
Specificity	This antibody detects endogenous levels of NIBAN at Human/Mouse/Rat	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	NIBAN	
Gene Name	FAM129A C1orf24 NIBAN GIG39	
Cellular localization	Cytoplasm . Membrane ; Lipid-anchor .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band		
Human Gene ID	116496	
Human Swiss-Prot Number	Q9BZQ8	
Alternative Names		
Background	This gene encodes a member of the family with sequence similarity 129 protein family. This gene is highly expressed in several cancer cells and may serve as a prognostic marker for certain cancers. The encoded protein may play a role in regulating p53-mediated apoptosis. [provided by RefSeq, Sep 2016],	




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkb

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C