

ASC rabbit pAb

Cat No.:ES1715

For research use only

Overview

Product Name	ASC rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300.
	Immunofluorescence: 1/200 - 1/1000. ELISA:
	1/40000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human ASC. AA range:10-59
Specificity	ASC Polyclonal Antibody detects endogenous levels
	of ASC protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	Apoptosis-associated speck-like protein containing a
	CARD
Gene Name	PYCARD
Cellular localization	Cytoplasm . Inflammasome . Endoplasmic
	reticulum . Mitochondrion . Nucleus . Upstream of
	caspase activation, a redistribution from the
	cytoplasm to the aggregates occurs. These appear as
	hollow, perinuclear spherical, ball-like structures
	(PubMed:11103777, PubMed:12191486,
	PubMed:15030775). Upon NLRP3 inflammasome
	activation redistributes to the perinuclear space
	localizing to endoplasmic reticulum and
	mitochondria (PubMed:12191486,
	PubMed:15030775). Localized primarily to the
	nucleus in resting monocytes/macrophages and
	rapidly redistributed to the cytoplasm upon
	pathogen infection (PubMed:19234215). Localized
	to large cytoplasmic aggregate appearing as a speck



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Purification	containing AIM2, PYCARD, CASP8 and bacterial DNA after infection with Francisella tularensis (By similarity); Golgi apparatus The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	21kD
Human Gene ID	29108
Human Swiss-Prot Number	Q9ULZ3
Alternative Names	PYCARD; ASC; CARD5; TMS1; Apoptosis-associated speck-like protein containing a CARD; hASC; Caspase recruitment domain-containing protein 5; PYD and CARD domain-containing protein; Target of methylation-induced silencing 1
Background	This gene encodes an adaptor protein that is composed of two protein-protein interaction domains: a N-terminal PYRIN-PAAD-DAPIN domain (PYD) and a C-terminal caspase-recruitment domain (CARD). The PYD and CARD domains are members of the six-helix bundle death domain-fold superfamily that mediates assembly of large signaling complexes in the inflammatory and apoptotic signaling pathways via the activation of caspase. In normal cells, this protein is localized to the cytoplasm; however, in cells undergoing apoptosis, it forms ball-like aggregates near the nuclear periphery. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C