

HRP–conjugated Goat Anti-Mouse IgG Catalog NO. : SA002 For research use only.

Overview		
Product name	HRP–conjugated Goat Anti-Mouse IgG	
Applications	WB,IHC,ELISA	
Recommended dilutions	WB IHC ELISA NOTE: Optimal dil	1:5000-1:10000 1:100-1:200 Use at an assay dependent concentration lutions should be determined by the end user.
Size	100µL	
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C for one year. Avoid repeated freeze- thaw cycles.	
Note	This product is for research use only,not for diagnostic or therapeutic use.	