

## CY3-conjugated Goat Anti-Mouse IgG Catalog NO.: SA007 For research use only.

Overview

Product name CY3-conjugated Goat Anti-Mouse IgG

Recommended dilutions 1:20-1:100 for most applications

NOTE: Optimal dilutions should be determined by the end user.

Size  $100\mu L$ 

**Storage** PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20°C for one year. Avoid repeated freeze-

thaw cycles.

Note This product is for research use only,not for diagnostic or

therapeutic use.

Fluorophore Cy3

Absorption Peak 550nm

Emission Peak 570nm