



ELK Biotechnology

MAKE RESEARCH EASIER

# Microsample ELISA kit

2000+ Kit

25µL only

New Arrivals



## Significant Sample Saving

**25µL** of sample required only! (traditional is 100µL), reduced by 75% directly. Precious clinical samples, model animal sera, and mouse sera with small blood collection volume will no longer be in short supply, allowing every sample of yours to be fully utilized!

## High Sensitivity

**pg/mL** Level LOD with high-affinity antibodies and an amplified system, enabling clear detection of even trace low-abundance targets. Uncover subtle biological changes and achieve more significant results for high-impact publications!

## High Specificity

Cross-reactivity **<5%**, so you can eliminate high background interference.

## Outstanding Precision

Intra- and inter-assay **CV < 8.5%**, delivering reliable and consistent results with excellent reproducibility. Trustworthy data provides a solid foundation for your outstanding research outcomes.

## ELK Biotechnology

### Research area of ELK Microsample ELISA kits

ELK Biotechnology Microsample ELISA Kits feature an extensive product line, successfully validated for over **800** hot-selling targets across **human, mouse, and rat species**, fully empowering your cross-disciplinary research.



- **Immunology and Inflammation Research**  
TNF-α, IL-6, IL-12, CRP, chemokines, IFN-γ, and other inflammatory mediators.
- **Cancer and Tumor**  
EGFR, CD19, CEA, PSA, and various growth factors, facilitating research in tumor mechanisms and diagnostics.
- **Metabolic Endocrinology**  
GLP1R, PPARα/δ, LEPR, PTH – enabling in-depth research on diabetes, obesity, and other metabolic diseases.
- **Cardiovascular Research**  
Cardiac Troponin (cTnI/T), BNP, PCSK9, NLRP3 – Key tools for precision diagnosis of cardiovascular diseases.
- **Stem Cell and Cell Therapy:** CD34, CD19, GPC3 – delivering critical data for regenerative medicine research.

Cat.NO	Product name	Detection range
ELK1157MS	Mouse IL6(Interleukin 6) Microsample ELISA Kit	7.82-500 pg/mL
ELK1156MS	Human IL6(Interleukin 6) Microsample ELISA Kit	7.82-500 pg/mL
ELK1387MS	Mouse TNFa(Tumor Necrosis Factor Alpha) Microsample ELISA Kit	15.63-1000 pg/mL
ELK1190MS	Human TNFa(Tumor Necrosis Factor Alpha) Microsample ELISA Kit	15.63-1000 pg/mL
ELK10900MS	Mouse TRAF6(TNF Receptor Associated Factor 6) Microsample ELISA Kit	0.32-20 µg/mL
ELK3775MS	Mouse sIL-6R(SolubleInterleukin 6 Receptor) Microsample ELISA Kit	0.32-20 ng/mL
ELK2507MS	Mouse IL-12p40(Interleukin 12 p40) Microsample ELISA Kit	15.63-1000 pg/mL
ELK9605MS	Mouse CXCL15(Chemokine C-X-C-Motif Ligand 15) Microsample ELISA Kit	78.13-5000 pg/mL
ELK1010MS	Mouse ECF(Eosinophil Chemotactic Factor) Microsample ELISA Kit	7.82-500 pg/mL
ELK10778MS	Mouse LECT2(Leukocyte Cell Derived Chemotaxin 2) Microsample ELISA Kit	1.57-100 ng/mL
ELK10579MS	Mouse SLC(Secondary Lymphoid Tissue Chemokine) Microsample ELISA Kit	31.25-2000 pg/mL
ELK2237MS	Mouse CD19(Cluster Of Differentiation 19) Microsample ELISA Kit	0.16-10 ng/mL
ELK2421MS	Mouse EGFR(Epidermal Growth Factor Receptor) Microsample ELISA Kit	93.75-6000 pg/mL
ELK9614MS	Mouse PSA(ProstateSpecific Antigen) Microsample ELISA Kit	0.16-10 ng/mL
ELK6119MS	Mouse GLP1R(Glucagon Like Peptide 1 Receptor) Microsample ELISA Kit	0.16-10 ng/mL
ELK2009MS	Mouse PPARa(Peroxisome Proliferator Activated Receptor Alpha) Microsample ELISA Kit	0.32-20 ng/mL
ELK9212MS	Rat hs-cTnI(High Sensitivity Cardiac Troponin I) Microsample ELISA Kit	15.63-1000 pg/mL
ELK6207MS	Mouse cTnT/TNNT2(Troponin T Type 2, Cardiac) Microsample ELISA Kit	15.63-1000 pg/mL
ELK2136MS	Mouse CTNnb1(Catenin Beta 1) Microsample ELISA Kit	0.16-10 ng/mL
ELK2237MS	Mouse CD19(Cluster Of Differentiation 19) Microsample ELISA Kit	0.16-10 ng/mL
ELK0918MS	Mouse GPC3(Glypican 3) Microsample ELISA Kit	0.16-10 ng/mL
ELK10393MS	Mouse LEPR(Leptin Receptor) Microsample ELISA Kit	3.13-200 ng/mL