


ELK Biotechnology


Transferrin Rabbit pAb Catalog NO.: EA188 For research use only.

Overview

Product name	Transferrin Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/5000-10000 Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Purified Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	77kDa
GenelD (Human)	7018
Human Swiss-Prot No.	P02787
Cellular localization	Secreted
Alternative Names	Apotransferrin antibody, beta 1 metal binding globulin antibody, serotransferrin antibody, TF antibody
Background	Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron in biological fluids. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization.



Western blot analysis of Human Blood sample with Transferrin Rabbit pAb EA188 diluted at 1:5,000.



Immunohistochemical analysis of paraffin-embedded Human Liver Tissue using Transferrin (EA188) Rabbit pAb diluted at 1:500.