



ELK Biotechnology


PI3 Kinase P85 β Rabbit pAb

Catalog NO.: EA210


For research use only.

Overview

Product name	PI3 Kinase P85 β Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human
Recommended dilutions	WesternBlot:1/1000-2000 Immunohistochemistry:1/50-100 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	82kDa
GeneID (Human)	5296
Human Swiss-Prot No.	O00459
Cellular localization	Cytosol,Nucleus
Alternative Names	P85 antibody, p85 beta antibody, PIK3R2 antibody, Phosphoinositide 3 kinase regulatory subunit 2 antibody
Background	PI3 Kinase P85 β mediating the association of the p110 catalytic unit to the plasma membrane.



Western blot analysis of HepG2 with PI3 Kinase P85β Rabbit pAb EA210 diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma using PI3 Kinase P85β (EA210) Rabbit pAb diluted at 1:200.