


ELK Biotechnology


c-Raf Rabbit pAb Catalog NO.: EA215 For research use only.

Overview

Product name	c-Raf Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human
Recommended dilutions	WesternBlot:1/1000-2000 Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	75kDa
GenelD (Human)	5894
Human Swiss-Prot No.	P04049
Cellular localization	N/A
Alternative Names	c-Raf proto-oncogene, serine/threonine kinase antibody, EC2.7.11.1 antibody, NS5 antibody, Raf 1 antibody
Background	A-Raf, B-Raf, and c-Raf (Raf-1) are the main effectors recruited by GTP- bound Ras to activate the MEK-MAP kinase pathway. C-Raf is a member of the Raf kinase family of serine/threonine-specific protein kinases, from the TKL (Tyrosine-kinase-like) group of kinases.



Western blot analysis of Hela with c-Raf Rabbit pAb EA215 diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma using c-Raf (EA215) Rabbit pAb diluted at 1:200.