



ELK Biotechnology


AMPK β 2 Rabbit pAb

Catalog NO.: EA231

For research use only.

Overview

Product name	AMPK β 2 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
GeneID (Human)	5565
Human Swiss-Prot No.	O43741
Cellular localization	Cytosol,Nucleus
Alternative Names	AMPK beta 2 antibody
Background	5' AMP-activated protein kinase or AMPK or 5' adenosine monophosphate-activated protein kinase is an enzyme that plays a role in cellular energy homeostasis.



Western blot analysis of 1) 293T, 2) HepG2, 3) Mouse Heart Tissue, 4) Rat Heart Tissue with AMPK β 2 Rabbit pAb EA231 diluted at 1:2,000.