


ELK Biotechnology Bombesin Receptor 2 Rabbit pAb Catalog NO.: EA263 For research use only.

Overview

Product name	Bombesin Receptor 2 Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50-70kDa
GenelD (Human)	2925
Human Swiss-Prot No.	P30550
Cellular localization	Cell membrane
Alternative Names	BB2 antibody
Background	Bombesin receptor 2 (BB2) is a member of a family of receptors that binds the 14 amino acid peptide, bombesin.



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with Bombesin Receptor 2 Rabbit pAb EA263 diluted at 1:2,000.

Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Bombesin Receptor 2 (EA263) Rabbit pAb diluted at 1:200.