



ELK Biotechnology


CNTF Rabbit pAb

Catalog NO.: EA279

For research use only.

Overview

Product name	CNTF Rabbit polyclonal antibody
Source	Rabbit
Applications	WB
Species reactivity	Human, Rat, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
GeneID (Human)	1270
Human Swiss-Prot No.	P26441
Cellular localization	Cytoplasm
Alternative Names	Ciliary neurotrophic factor, HCNTF
Background	CNTF is a survival factor for various neuronal cell types. Seems to prevent the degeneration of motor axons after axotomy.



Western blot analysis of 1) Hela, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with CNTF Rabbit pAb EA279 diluted at 1:2,000.