



ELK Biotechnology


GRIK2(GluR6) Rabbit pAb

Catalog NO.: EA297

For research use only.

Overview

Product name	GRIK2(GluR6) Rabbit polyclonal antibody
Source	Rabbit
Applications	IHC
Species reactivity	Human
Recommended dilutions	Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Synthetic Peptide
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	103kDa
GenelD (Human)	2898
Human Swiss-Prot No.	Q13002
Cellular localization	Cell junction, Cell membrane, Membrane, Postsynaptic cell membrane, Synapse
Alternative Names	Excitatory amino acid receptor 4, GLR6, GLUK2, Glutamate receptor 6
Background	Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system.



Immunohistochemical analysis of paraffin-embedded Human Colon Tissue using GRIK2(GluR6) (EA297) Rabbit pAb diluted at 1:200.