


## **ELK Biotechnology**

c-Jun Rabbit pAb Catalog NO.: EA329 For research use only.

Overview

Product name	c-Jun Rabbit polyclonal antibody
Source	Rabbit
Applications	WB, IHC
Species reactivity	Human, Mouse
Recommended dilutions	WesternBlot:1/1000-2000 Immunohistochemistry:1/100-200 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	lgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	43kDa
GenelD (Human)	3725
Human Swiss-Prot No.	P05412
Cellular localization	Nucleus
Alternative Names	Activator Protein 1 antibody, AP1 antibody, p39 antibody
Background	c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1).



14KD

Immunohistochemical analysis of paraffin-embedded Human Ovarain Carcinoma Tissue using c-Jun (EA329) Rabbit pAb diluted at 1:200.

Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using c-Jun (EA329) Rabbit pAb diluted at 1:200

Western blot analysis of Mouse Brain Tissue Lysate using c-Jun (EA329) Rabbit pAb diluted at 1:2000.