



ELK Biotechnology


Gab1 Rabbit pAb

Catalog NO.: EA360

For research use only.

Overview

Product name	Gab1 Rabbit polyclonal antibody
Source	Rabbit
Applications	IHC
Species reactivity	Human, Mouse, Rat
Recommended dilutions	Immunohistochemistry:1/200-500 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	Human
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	76kDa
GeneID (Human)	2549
Human Swiss-Prot No.	Q13480
Cellular localization	Cytosol
Alternative Names	GRB 2 associated binder 1 antibody
Background	The Gab1 protein is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Gab1 (EA360) Rabbit pAb diluted at 1:500