

PKHG1 rabbit pAb

Cat No.: ES10011

For research use only

Overview

Product Name PKHG1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 780-860

Specificity PKHG1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Pleckstrin homology domain-containing family G

member 1

Gene Name PLEKHG1 KIAA1209

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 152kD
Human Gene ID 57480
Human Swiss-Prot Number Q9ULL1

Alternative Names

Background PTM:Phosphorylated upon DNA damage, probably

by ATM or ATR., similarity: Contains 1 DH

(DBL-homology) domain., similarity: Contains 1 PH

domain.,