

KCNS1 rabbit pAb

Cat No.: ES10033

For research use only

Overview

Product Name KCNS1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 470-550

Specificity KCNS1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Potassium voltage-gated channel subfamily S
member 1 (Delayed-rectifier K(+) channel alpha

subunit 1) (Voltage-gated potassium channel

subunit Kv9.1)

Gene Name KCNS1

Cellular localization Cell membrane; Multi-pass membrane protein.

May not reach the plasma membrane but remain in an intracellular compartment in the absence of

KCNB1 or KCNB2 (PubMed:10484328). .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.


Clonality Polyclonal
Concentration 1 mg/ml
Observed band 57kD
Human Gene ID 3787
Human Swiss-Prot Number O96KK3

Alternative Names

Background Voltage-gated potassium channels form the largest

and most diversified class of ion channels and are present in both excitable and nonexcitable cells.





Their main functions are associated with the regulation of the resting membrane potential and the control of the shape and frequency of action potentials. The alpha subunits are of 2 types: those that are functional by themselves and those that are electrically silent but capable of modulating the activity of specific functional alpha subunits. The protein encoded by this gene is not functional by itself but can form heteromultimers with member 1 and with member 2 (and possibly other members) of the Shab-related subfamily of potassium voltage-gated channel proteins. This gene belongs to the S subfamily of the potassium channel family. [provided by RefSeq, Jul 2008],

Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000