

## **POTEE** rabbit pAb

Cat No.: ES10052

For research use only

## Overview

Product Name POTEE rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 350-430

**Specificity** POTEE Polyclonal Antibody detects endogenous

levels of protein.

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at  $-20^{\circ}$ C. Avoid repeated freeze-thaw cycles.

**Protein Name** POTE ankyrin domain family member E

(ANKRD26-like family C member 1A) (Prostate, ovary, testis-expressed protein on chromosome 2)

(POTE-2)

Gene Name POTEE A26C1A POTE2

**Cellular localization** extracellular space, extracellular exosome, blood

microparticle,

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 118kD
Human Gene ID 445582
Human Swiss-Prot Number Q6S8J3

**Alternative Names** 

**Background** similarity:Belongs to the actin

family., similarity: Contains 5 ANK

repeats., similarity: In the C-terminal section; belongs

to the actin family., similarity: In the N-terminal

section; belongs to the POTE family.,