

PPM1H rabbit pAb

Cat No.: ES10058

For research use only

Overview

Product Name PPM1H rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 290-370

Specificity PPM1H Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Protein phosphatase 1H (EC 3.1.3.16)

Gene Name PPM1H ARHCL1 KIAA1157 URCC2

Cellular localization Nucleus . Cytoplasm .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 56kD
Human Gene ID 57460
Human Swiss-Prot Number Q9ULR3

Alternative Names

Background catalytic activity: A phosphoprotein + H(2)O = a

protein + phosphate., similarity: Belongs to the PP2C family., similarity: Contains 1 PP2C-like domain.,