

Cleaved-Caspase-9 p35 (D315) rabbit pAb

Cat No.: ES1009

For research use only

Overview

Product Name Cleaved-Caspase-9 p35 (D315) rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA **Species Cross-Reactivity** Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300

Immunogen The antiserum was produced against synthesized

range:266-315

Specificity Cleaved-Caspase-9 p35 (D315) Polyclonal Antibody

detects endogenous levels of fragment of activated Caspase-9 p35 protein resulting from cleavage

peptide derived from human Caspase 9. AA

adjacent to D315.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Caspase9
Gene Name CASP9

Cellular localizationnucleus,mitochondrion,cytosol,apoptosome,PurificationThe antibody was affinity-purified from rabbit
antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band35 46kDHuman Gene ID842Human Swiss-Prot NumberP55211

Alternative Names CASP9; MCH6; Caspase-9; CASP-9; Apoptotic

protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6;

ICE-LAP6

Background CASP9 encodes a member of the cysteine-aspartic

acid protease (caspase) family. Sequential activation

of caspases plays a central role in the



execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. Caspase 9 can undergo autoproteolytic processing and activation by the apoptosome, a protein complex of cytochrome c and the apoptotic peptidase activating factor 1; this step is thought to be one of the earliest in the caspase activation cascade. Caspase 9 is thought to play a central role in apoptosis and to be a tumor suppressor. Alternative splicing results in multiple transcript variants.