

## RB40A rabbit pAb

Cat No.: ES10127

For research use only

## Overview

Product Name RB40A rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

**Immunogen** Synthesized peptide derived from human protein .

at AA range: 160-240

**Specificity** RB40A Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ras-related protein Rab-40A (SOCS box-containing

protein RAR2A) (Protein Rar-2)

Gene Name RAB40A

Cellular localizationCell membrane ; Lipid-anchor ; Cytoplasmic side .PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.


Clonality Polyclonal
Concentration 1 mg/ml
Observed band 30kD
Human Gene ID 142684
Human Swiss-Prot Number Q8WXH6

**Alternative Names** 

**Background** This gene encodes a member of the Rab40

subfamily of Rab small GTP-binding proteins that contains a C-terminal suppressors of cytokine signaling box. [provided by RefSeq, Apr 2010],





Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night