

## **RAB9A** rabbit pAb


## Cat No.:ES10135

For research use only

## Overview

Product Name Host species Applications Species Cross-Reactivity Recommended dilutions Immunogen	RAB9A rabbit pAb Rabbit WB;ELISA Human;Rat;Mouse WB 1:500-2000 ELISA 1:5000-20000 Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	RAB9A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ras-related protein Rab-9A
Gene Name	RAB9A RAB9
Cellular localization	Cell membrane ; Lipid-anchor ; Cytoplasmic side . Endoplasmic reticulum membrane . Golgi apparatus membrane . Late endosome . Cytoplasmic vesicle, phagosome membrane ; Lipid-anchor ; Cytoplasmic side . Cytoplasmic vesicle, phagosome . Cytoplasmic vesicle
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	22kD
Human Gene ID	9367
Human Swiss-Prot Number	P51151
Alternative Names	
Background	function:Involved in the transport of proteins between the endosomes and the trans Golgi network.,similarity:Belongs to the small GTPase superfamily. Rab family.,





Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night