

## RL1D1 rabbit pAb

Cat No.: ES10180

For research use only

## Overview

Product Name RL1D1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 280-360

**Specificity** RL1D1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20℃. Avoid repeated freeze-thaw cycles.Protein NameRibosomal L1 domain-containing protein 1 (CATX-11)

(Cellular senescence-inhibited gene protein) (Protein

PBK1)

Gene Name RSL1D1 CATX11 CSIG PBK1 L12

**Cellular localization** Nucleus, nucleolus . Colocalizes with ING1 in the

nucleolus after UV stress. .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band53kDHuman Gene ID26156Human Swiss-Prot Number076021

**Alternative Names** 

Background similarity: Belongs to the ribosomal protein L1P

family. Highly divergent., tissue specificity: Placenta.,