

ELL3 rabbit pAb

Cat No.: ES10187

For research use only

Overview

Product Name ELL3 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 210-290

Specificity ELL3 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name RNA polymerase II elongation factor ELL3

Gene Name ELL3
Cellular localization Nucleus .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 43kD
Human Gene ID 80237
Human Swiss-Prot Number Q9HB65


Alternative Names

Background function: Elongation factor that can increase the

catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.,similarity:Belongs to the ELL/occludin family.,tissue specificity:Testis

specific.,





Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night