

RPAP3 rabbit pAb

Cat No.: ES10189

For research use only

Overview

Product Name RPAP3 rabbit pAb

Host species Rabbit **Applications** WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 280-360

RPAP3 Polyclonal Antibody detects endogenous Specificity

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20° C. Avoid repeated freeze-thaw cycles. **Storage**

RNA polymerase II-associated protein 3 **Protein Name**

RPAP3 **Gene Name**

Cellular localization R2TP complex,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml **Observed band** 73kD **Human Gene ID** 79657 **Human Swiss-Prot Number** Q9H6T3

Alternative Names

Background This gene encodes an RNA polymerase II-associated

> protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar

2009],