

SCN1B rabbit pAb

Cat No.: ES10279

For research use only

Overview

Product Name SCN1B rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity SCN1B Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Sodium channel subunit beta-1

Gene Name SCN1B

Cellular localization [Isoform 1]: Cell membrane; Single-pass type I

membrane protein . Perikaryon . Cell projection . Cell projection, axon . Detected at nodes of Ranvier on the sciatic nerve. .; [Isoform 2]: Perikaryon . Cell projection . Secreted . Detected on Purkinje cel The antibody was affinity-purified from rabbit

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 23kD
Human Gene ID 6324
Human Swiss-Prot Number 007699

Alternative Names


Background Voltage-gated sodium channels are heteromeric

proteins that function in the generation and propagation of action potentials in muscle and neuronal cells. They are composed of one alpha and

two beta subunits, where the alpha subunit







provides channel activity and the beta-1 subunit modulates the kinetics of channel inactivation. This gene encodes a sodium channel beta-1 subunit. Mutations in this gene result in generalized epilepsy with febrile seizures plus, Brugada syndrome 5, and defects in cardiac conduction. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2009],

Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night



23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C

+86-27-59760950