

## TLX3 rabbit pAb

Cat No.: ES10347

For research use only

## Overview

Product Name TLX3 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

**Immunogen** Synthesized peptide derived from part region of

human protein

**Specificity** TLX3 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name T-cell leukemia homeobox protein 3 (Homeobox

protein Hox-11L2)

Gene Name TLX3 HOX11L2

Cellular localization Nucleus .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 32kD
Human Gene ID 30012
Human Swiss-Prot Number O43711

**Alternative Names** 

**Background** The protein encoded by this gene is an orphan

homeobox protein that encodes a DNA-binding nuclear transcription factor. A translocation

[t(5;14)(q35;q32)] involving this gene is associated with T-cell acute lymphoblastic leukemia (T-ALL) in children and young adults. [provided by RefSeq, Nov


2015],



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

