

TENN rabbit pAb

Cat No.: ES10350

For research use only

Overview

Product Name TENN rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity TENN Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Tenascin-N (TN-N)

Gene Name TNN

Cellular localization Secreted, extracellular space, extracellular matrix .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 142kD
Human Gene ID 63923
Human Swiss-Prot Number Q9UQP3

Alternative Names

Background function:Involved in neurite outgrowth and cell

migration in hippocampal

explants., similarity: Belongs to the tenascin

family., similarity: Contains 1 fibrinogen C-terminal

domain., similarity: Contains 3 EGF-like

domains., similarity: Contains 9 fibronectin type-III

domains.,



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com