

UBP3 rabbit pAb

Cat No.: ES10416

For research use only

Overview

Product Name UBP3 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity UBP3 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Ubiquitin carboxyl-terminal hydrolase 3 (EC

3.4.19.12) (Deubiquitinating enzyme 3) (Ubiquitin thioesterase 3) (Ubiquitin-specific-processing

protease 3)

Gene Name USP3

Cellular localization Nucleus . Localizes preferentially with

monoubiquitinated H2A to chromatin.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 57kD
Human Gene ID 9960
Human Swiss-Prot Number Q9Y6I4

Alternative Names

Background catalytic activity: Ubiquitin C-terminal thioester +

H(2)O = ubiquitin + a thiol., similarity: Belongs to the peptidase C19 family., similarity: Contains 1 UBP-type zinc finger., tissue specificity: Expressed in all tissues examined, with strongest expression in pancreas.,



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com