

## **UBP31** rabbit pAb

Cat No.: ES10417

For research use only

## Overview

Product Name UBP31 rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity UBP31 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ubiquitin carboxyl-terminal hydrolase 31 (EC

3.4.19.12) (Deubiquitinating enzyme 31) (Ubiquitin thioesterase 31) (Ubiquitin-specific-processing

protease 31)

Gene Name USP31 KIAA1203

Cellular localization nucleus,

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band148kDHuman Gene ID57478Human Swiss-Prot NumberQ70CQ4

**Alternative Names** 

Background catalytic activity: Ubiquitin C-terminal thioester +

H(2)O = ubiquitin + a thiol.,function:May recognize and hydrolyze the peptide bond at the C-terminal Gly of ubiquitin. Involved in the processing of poly-ubiquitin precursors as well as that of ubiquinated proteins.,similarity:Belongs to the



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



peptidase C19 family., tissue specificity: Widely expressed.,



+86-27-59760950