

UBE20 rabbit pAb

Cat No.: ES10440

For research use only

Overview

Product Name UBE2O rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity UBE2O Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ubiquitin-conjugating enzyme E2 O (EC 6.3.2.19)

(Ubiquitin corrier protein O) (Ubiquitin conjugating

(Ubiquitin carrier protein O) (Ubiquitin-conjugating enzyme E2 of 230 kDa) (Ubiquitin-conjugating enzyme E2-230K) (Ubiquitin-protein ligase O)

Gene Name UBE20 KIAA1734

Cellular localization Cytoplasm . Nucleus . Mainly localizes to the

cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 142kD
Human Gene ID 63893
Human Swiss-Prot Number Q9C0C9

Alternative Names

Background catalytic activity:ATP + ubiquitin + protein lysine =

AMP + diphosphate + protein

N-ubiquityllysine., function: Catalyzes the covalent

attachment of ubiquitin to other

proteins.,pathway:Protein modification; protein



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



ubiquitination., similarity: Belongs to the ubiquitin-conjugating enzyme family., tissue specificity: Predominantly expressed in skeletal muscle and heart.,



+86-27-59760950