

## **UB2Q2** rabbit pAb

Cat No.: ES10442

For research use only

## Overview

Product Name UB2Q2 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

**Immunogen** Synthesized peptide derived from part region of

human protein

**Specificity** UB2Q2 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Ubiquitin-conjugating enzyme E2 Q2 (EC 6.3.2.19)

(Ubiquitin carrier protein Q2) (Ubiquitin-protein

(Ubiquitin carrier protein Q2) (Ubiquitin-protein

ligase Q2)

**Gene Name** UBE2Q2 **Cellular localization** Cytoplasm .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 41kD
Human Gene ID 92912
Human Swiss-Prot Number Q8WVN8

**Alternative Names** 

Background catalytic activity:ATP + ubiquitin + protein lysine =

AMP + diphosphate + protein

N-ubiquityllysine.,function:Catalyzes the covalent


attachment of ubiquitin to other

proteins., pathway: Protein modification; protein

ubiquitination., similarity: Belongs to the ubiquitin-conjugating enzyme family., tissue







specificity:Detected in hypopharyngeal head and neck squamous cell carcinoma, in tumor masses and invasive epithelium.,

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950