

## Cleaved-MMP-27 (Y99) rabbit pAb

Cat No.: ES1048

For research use only

## Overview

Cleaved-MMP-27 (Y99) rabbit pAb **Product Name** 

**Host species** Rabbit WB;ELISA **Applications** 

**Species Cross-Reactivity** Human; Monkey

**Recommended dilutions** Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not

yet tested in other applications.

The antiserum was produced against synthesized **Immunogen** 

peptide derived from human MMP27. AA

range:80-129

**Specificity** Cleaved-MMP-27 (Y99) Polyclonal Antibody detects

> endogenous levels of fragment of activated MMP-27 protein resulting from cleavage adjacent to Y99.

**Formulation** Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles. Storage

**Protein Name** Matrix metalloproteinase-27

**Gene Name** MMP27

Cellular localization Endoplasmic reticulum membrane; Peripheral

membrane protein . Retained in the endoplasmic

reticulum...

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal 1 mg/ml Concentration **Observed band** 48kD **Human Gene ID** 64066 **Human Swiss-Prot Number** Q9H306

**Alternative Names** 

MMP27; Matrix metalloproteinase-27; MMP-27 Background Proteins of the matrix metalloproteinase (MMP)


family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue




+86-27-59760950 ELKbio@ELKbiotech.com



remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. [provided by RefSeq, Jul 2008],





MMP27 -- 48
-- 34
-- 26
-- 19

(kD)

Western blot analysis of lysates from COS7 cells, treated with etoposide 25uM 1h, using MMP27 (Cleaved-Tyr99) Antibody. The lane on the right is blocked with the synthesized peptide.

