

ZBTB2 rabbit pAb

Cat No.:ES10480

For research use only

Overview

Product Name	ZBTB2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
	human protein
Specificity	ZBTB2 Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger and BTB domain-containing protein 2
Gene Name	ZBTB2 KIAA1483 ZNF437
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	56kD
Human Gene ID	57621
Human Swiss-Prot Number	Q8N680
Alternative Names	
Background	function:May be involved in transcriptional
	regulation.,PTM:Phosphorylated upon DNA damage,
	probably by ATM or ATR., similarity: Contains 1 BTB
	(POZ) domain., similarity: Contains 4 C2H2-type zinc
	fingers.,



ELKbio@ELKbiotech.com

m www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C