

T106C rabbit pAb

Cat No.:ES10641

For research use only

Overview

Product Name	T106C rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	T106C Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Transmembrane protein 106C (Endoplasmic
	reticulum membrane protein overexpressed in
	cancer)
Gene Name	TMEM106C EMOC
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass
	membrane protein . Membrane ; Lipid-anchor .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	27kD
Human Gene ID	79022
Human Swiss-Prot Number	Q9BVX2
Alternative Names	
Background	similarity:Belongs to the TMEM106 family.,




ELKbio@ELKbiotech.com

om www.elk

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C