

BT3A1 rabbit pAb

Cat No.: ES10865

For research use only

Overview

Product Name BT3A1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity BT3A1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles. Protein Name Butyrophilin subfamily 3 member A1 (CD antigen

CD277)

Gene Name BTN3A1 BTF5

Cell ular localization Cell membrane ; Single-pass type I membrane

protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 56kD
Human Gene ID 11119
Human Swiss-Prot Number 000481

Alternative Names

Background The butyrophilin (BTN) genes are a group of major

histocompatibility complex (MHC)-associated genes that encode type I membrane proteins with 2 extracellular immunoglobulin (Ig) domains and an intracellular B30.2 (PRYSPRY) domain. Three

subfamilies of human BTN genes are located in the MHC class I region: the single-copy BTN1A1 gene



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



(MIM 601610) and the BTN2 (e.g., BTN2A1; MIM 613590) and BTN3 (e.g., BNT3A1) genes, which have undergone tandem duplication, resulting in 3 copies of each (summary by Smith et al., 2010 [PubMed 20208008]).[supplied by OMIM, Nov 2010],



+86-27-59760950