

SKI rabbit pAb

Cat No.: ES10887

For research use only

Overview

Product Name SKI rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity SKI Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Ski oncogene (Proto-oncogene c-Ski)

Gene Name SKI

Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 80kD
Human Gene ID 6497
Human Swiss-Prot Number P12755

Alternative Names

Background SKI proto-oncogene(SKI) Homo sapiens This

gene encodes the nuclear protooncogene protein homolog of avian sarcoma viral (v-ski) oncogene. It functions as a repressor of TGF-beta signaling, and may play a role in neural tube development and muscle differentiation. [provided by RefSeq, Oct

2009],



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com