

SP3 rabbit pAb

Cat No.: ES10910

For research use only

Overview

Product Name SP3 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity SP3 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Transcription factor Sp3 (SPR-2)

Gene Name SP3

Cellular localization Nucleus. Nucleus, PML body. Localizes to the nuclear

periphery and in nuclear dots when sumoylated.

Some localization in PML nuclear bodies.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 85kD
Human Gene ID 6670
Human Swiss-Prot Number Q02447

Alternative Names

Background This gene belongs to a family of Sp1 related genes

that encode transcription factors that regulate transcription by binding to consensus GC- and GT-box regulatory elements in target genes. This protein contains a zinc finger DNA-binding domain and several transactivation domains, and has been reported to function as a bifunctional transcription



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



factor that either stimulates or represses the transcription of numerous genes. Transcript variants encoding different isoforms have been described for this gene, and one has been reported to initiate translation from a non-AUG (AUA) start codon. Additional isoforms, resulting from the use of alternate downstream translation initiation sites, have also been noted. A related pseudogene has been identified on chromosome 13. [provided by RefSeq, Feb 2010],



+86-27-59760950