

FAM3C rabbit pAb

Cat No.:ES10951

For research use only

Overview

Product Name	FAM3C rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human;Mouse;Rat	
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000	
Immunogen	Synthesized peptide derived from part region of	
	human protein	
Specificity	FAM3C Polyclonal Antibody detects endogenous	
	levels of protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Protein FAM3C (Interleukin-like EMT inducer)	
Gene Name	FAM3C ILEI GS3786	
Cellular localization	Secreted . Cytoplasmic vesicle . Cytoplasmic in some	
	cancer cells.	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	24kD	
Human Gene ID	10447	
Human Swiss-Prot Number	Q92520	
Alternative Names		
Background	This gene is a member of the family with sequence	
	similarity 3 (FAM3) family and encodes a secreted	
	protein with a GG domain. A change in expression of	
	this protein has been noted in pancreatic	
	cancer-derived cells. [provided by RefSeq, Mar	
	2010],	




ELKbio@ELKbiotech.com

n www.elkb

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C