

REL1 rabbit pAb

Cat No.: ES10982

For research use only

Overview

Product Name REL1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity REL1 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Prorelaxin H1 [Cleaved into: Relaxin B chain; Relaxin

A chain]

Gene Name RLN1
Cellular localization Secreted.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 20kD
Human Gene ID 6013
Human Swiss-Prot Number P04808

Alternative Names

Background relaxin 1(RLN1) Homo sapiens Relaxins are

known endocrine and autocrine/paracrine hormones, belonging to the insulin gene


superfamily. In humans there are three non-allelic relaxin genes, RLN1, RLN2 and RLN3, where RLN1 and RLN2 share high sequence homology. The protein encoded by this gene is synthesized as a single-chain polypeptide but the active form consists



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





of an A chain and a B chain linked by disulfide bonds. Relaxin is produced by the ovary, and targets the mammalian reproductive system to ripen the cervix, elongate the pubic symphysis and inhibit uterine contraction. It may have additional roles in enhancing sperm motility, regulating blood pressure, controlling heart rate and releasing oxytocin and vasopressin. [provided by RefSeq, Jan 2013],

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

