

ADA22 rabbit pAb

Cat No.: ES11016

For research use only

Overview

Product Name ADA22 rabbit pAb

Host species Rabbit
Applications WB;ELISA
Species Cross-Reactivity Human;Mouse

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity ADA22 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Disintegrin and metalloproteinase

domain-containing protein 22 (ADAM 22) (Metalloproteinase-disintegrin ADAM22-3) (Metalloproteinase-like, disintegrin-like, and

cysteine-rich protein 2)

Gene Name ADAM22 MDC2

Cell ular localization Cell membrane ; Single-pass type I membrane

protein. Cell projection, axon.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 99kD
Human Gene ID 53616
Human Swiss-Prot Number Q9P0K1

Alternative Names

Background ADAM metallopeptidase domain 22(ADAM22) Homo

sapiens This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain)

family. Members of this family are



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. Unlike other members of the ADAM protein family, the protein encoded by this gene lacks metalloprotease activity since it has no zinc-binding motif. This gene is highly expressed in the brain and may function as an integrin ligand in the brain. In mice, it has been shown to be essential for correct myelination in the peripheral nervous system. Alternative splicing results in several transcript variants. [provided by RefSeq, Dec 2010],



+86-27-59760950