

## ATL1 rabbit pAb

## Cat No.:ES11143

For research use only

## Overview

Product Name	ATL1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
	human protein
Specificity	ATL1 Polyclonal Antibody detects endogenous levels
	of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$ . Avoid repeated freeze-thaw cycles.
Protein Name	ADAMTS-like protein 1 (ADAMTSL-1) (Punctin-1)
Gene Name	ADAMTSL1 ADAMTSR1 C9orf94 UNQ528/PRO1071
Cellular localization	Secreted, extracellular space, extracellular matrix .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	193kD
Human Gene ID	92949
Human Swiss-Prot Number	Q8N6G6
Alternative Names	
Background	ADAMTS like 1(ADAMTSL1) Homo sapiens This
	gene encodes a secreted protein and member of the
	ADAMTS (a disintegrin and metalloproteinase with
	thrombospondin motif) family. This protein lacks the
	metalloproteinase and disintegrin-like domains,
	which are typical of the ADAMTS family, but contains
	other ADAMTS domains, including the
	thrombospondin type 1 motif. This protein may have
	important functions in the extracellular matrix.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Jul 2008],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C