

## TIMD4 rabbit pAb

## Cat No.:ES11159

For research use only

## Overview

Product Name	TIMD4 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of
	human protein
Specificity	TIMD4 Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	T-cell immunoglobulin and mucin
	domain-containing protein 4 (TIMD-4) (T-cell
	membrane protein 4) (TIM-4)
Gene Name	TIMD4 TIM4
Cellular localization	Cell membrane ; Single-pass type I membrane
	protein . Secreted, extracellular exosome .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	41kD
Human Gene ID	91937
Human Swiss-Prot Number	Q96H15
Alternative Names	
Background	caution:The sequence shown here is derived from
	an Ensembl automatic analysis pipeline and should
	be considered as preliminary
	data.,domain:Recognizes phosphatidyl serine via its
	immunoglobulin
	domain.,function:Phosphatidylserine receptor that



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



enhances the engulfment of apoptotic cells. Involved in regulating T-cell proliferation and lymphotoxin signaling. Ligand for HAVCR1/TIMD1.,similarity:Belongs to the immunoglobulin superfamily. TIM family.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C