

NAAA rabbit pAb

Cat No.:ES11250

For research use only

Overview

Product Name	NAAA rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein .
	at AA range: 230-310
Specificity	NAAA Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	N-acylethanolamine-hydrolyzing acid amidase (EC
	3.5.1) (Acid ceramidase-like protein)
	(N-acylsphingosine amidohydrolase-like) (ASAH-like
	protein) [Cleaved into:
	N-acylethanolamine-hydrolyzing acid a
Gene Name	NAAA ASAHL PLT
Cellular localization	Lysosome . Membrane ; Peripheral membrane
	protein .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	39kD
Human Gene ID	27163
Human Swiss-Prot Number	Q02083
Alternative Names	
Background	This gene encodes an
	N-acylethanolamine-hydrolyzing enzyme which is
	highly similar to acid ceramidase. Multiple transcript
	variants encoding different isoforms have been




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





found for this gene. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C