

GLMN rabbit pAb

Cat No.:ES11319

For research use only

Overview

Product Name	GLMN rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at
	AA range: 80-160
Specificity	GLMN Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	Glomulin (FK506-binding protein-associated protein)
	(FAP) (FKBP-associated protein)
Gene Name	GLMN FAP48 FAP68 VMGLOM
Cellular localization	intracellular,cytoplasm,cullin-RING ubiquitin ligase
	complex,Cul2-RING ubiquitin ligase
	complex,Cul3-RING ubiquitin ligase
	complex,Cul4A-RING E3 ubiquitin ligase complex,
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	65kD
Human Gene ID	11146
Human Swiss-Prot Number	Q92990
Alternative Names	
Background	glomulin, FKBP associated protein(GLMN) Homo
	sapiens This gene encodes a phosphorylated
	protein that is a member of a Skp1-Cullin-F-box-like
	complex. The protein is essential for normal
	development of the vasculature and mutations in



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



this gene have been associated with glomuvenous malformations, also called glomangiomas. Multiple splice variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C