

## **BACE2** rabbit pAb

## Cat No.:ES11348

For research use only

## Overview

Product Name	BACE2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human; Rat; Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at
-	AA range: 210-290
Specificity	BACE2 Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!{ m C}$ . Avoid repeated freeze-thaw cycles.
Protein Name	Beta-secretase 2 (EC 3.4.23.45) (Aspartic-like
	protease 56 kDa) (Aspartyl protease 1) (ASP1) (Asp
	1) (Beta-site amyloid precursor protein cleaving
	enzyme 2) (Beta-site APP cleaving enzyme 2) (Down
	reg
Gene Name	BACE2 AEPLC ALP56 ASP21 CDA13 UNQ418/PRO852
Cellular localization	Cell membrane ; Single-pass type I membrane
	protein . Golgi apparatus . Endoplasmic reticulum.
	Endosome.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	56kD
Human Gene ID	25825
Human Swiss-Prot Number	Q9Y5Z0
Alternative Names	
Background	This gene encodes an integral membrane
	glycoprotein that functions as an aspartic protease.
	The encoded protein cleaves amyloid precursor



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



protein into amyloid beta peptide, which is a critical step in the etiology of Alzheimer's disease and Down syndrome. The protein precursor is further processed into an active mature peptide. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C