



PLAP rabbit pAb

Cat No.:ES11397

For research use only

Overview

Product Name	PLAP rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 531-580
Specificity	PLAP Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Phospholipase A-2-activating protein (PLA2P) (PLAP)
Gene Name	PLAA PLAP
Cellular localization	Nucleus . Cytoplasm . Cell junction, synapse . Recruited to damaged lysosomes decorated with K48-linked ubiquitin chains. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	87kD
Human Gene ID	9373
Human Swiss-Prot Number	Q9Y263
Alternative Names	
Background	function:Plays an important role in the regulation of specific inflammatory disease processes.,sequence caution:Translated as Gln.,similarity:Belongs to the WD repeat PLAP family.,similarity:Contains 1 PFU domain.,similarity:Contains 1 PUL domain.,similarity:Contains 7 WD repeats.,

