

ZBP1 rabbit pAb

Cat No.: ES11422

For research use only

Overview

Product Name ZBP1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 21-70

Specificity ZBP1 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Z-DNA-binding protein 1 (Tumor stroma and

activated macrophage protein DLM-1)

Gene Name ZBP1 C20orf183 DLM1

Cellular localization Cytoplasm . Nucleus . Mainly cytoplasmic

(PubMed:16876127, PubMed:16990255).

Accumulates in the nucleus in response to influenza A virus (IAV) infection: senses IAV defective viral genomes RNA in the nucleus (By similarity). .;

[Isoform 7]: Cytoplasm . Nu

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 47kD
Human Gene ID 81030
Human Swiss-Prot Number O9H171

Alternative Names

Background This gene encodes a Z-DNA binding protein. The


encoded protein plays a role in the innate immune response by binding to foreign DNA and inducing



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





type-I interferon production. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011],

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950