

SEM4C rabbit pAb

Cat No.: ES11818

For research use only

Overview

Product Name SEM4C rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity SEM4C Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Semaphorin-4C

Gene Name SEMA4C KIAA1739 SEMAI UNQ5855/PRO34487
Cellular localization Cell junction, synapse, postsynaptic density

membrane; Single-pass type I membrane protein. Cytoplasmic vesicle, secretory vesicle, synaptic vesicle membrane; Single-pass type I membrane

protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 91kD
Human Gene ID 54910
Human Swiss-Prot Number Q9C0C4

Alternative Names

Background similarity:Belongs to the semaphorin

family., similarity: Contains 1 Ig-like C2-type

(immunoglobulin-like) domain., similarity: Contains 1

PSI domain., similarity: Contains 1 Sema

domain., subunit: Interacts with GIPC PDZ domain.



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



Interacts with NCDN.,

